

our ongoing work [3] directed to manage together the diversity and flexibility within a pool of ligand candidates for bioassays, [4] we present here the results of a study concerning the replacement of CH₂ group (X in figure) in Irurre compounds, [2] with oxygen, sulfur, or with a SO₂ group. The structural properties of the resulting molecules were studied in the solid state, by single crystal X-ray diffraction, and calculated in the gas phase, by ab-initio methods. In each case the energy barrier to be overcome for the enantiomers interconversion as well as the transition state have been determined. The resulting scale of flexibility has been correlated with the chemical and structural features of the diverse library members.

[1] Evans, B.E., et al., *J. Med. Chem.*, 1988, 31, 2235.

[2] Irurre J., et al., *Can. J. Chem.*, 1994, 72, 334.

[3] Altamura, M., et al., *Tetrahedron* 2006, 62, 6754 and references herein.

[4] Guidi, A., WO 2006097449. Chem Abstr. 2006, 145, 356813.

MS20 P10

Quantification of pharmaceuticals in solid dosage forms Hana Petrickova^a ^aZentiva a.s., R&D Analytical department, U kabelovny 130, 10237 Prague 10 – Dolni Mecholupy, Czech Republic.

E-mail: hana.petrickova@zentiva.cz

Keywords: quantitative XRPD, pharmaceuticals, phase identification

XRPD is nowadays a routine widespread tool for characterisation pharmaceutical solids. The qualitative phase analysis is essential for development either API itself or the final solid dosage form and a lot of applications were introduced ranging polymorphic screenings, pre-formulation, formulation, stability testing or crystallography [1]. On the other hand importance of quantitative phase analysis (QPA) of polycrystalline mixtures comes into ever-increasing attention. After identification comes, of course, question: "How much?". Crucial for the analysis is to decide which of the quantitative information is expected, concrete number, limit test or simple positive/negative evidence.

Examples of different approaches of QPA will be presented using the FullPat [2], the Rietveld [3] and the single peak method. Influence of the sample (amorphous/crystalline, with/ without knowledge of 3D crystal structure, crystal size and shape) will be discussed on two examples: crystalline three component mixture and almost amorphous tableting mixture.

[1] A. and M. Zakrzewski: Solid State Charact. of Pharmaceuticals. Pergamon, Danbury, Connecticut, 2006.

[2] Chipera S.J., Bish D.L., *J. Appl. Cryst.*, 2002, 35, 744.

[3] Iyengar S.S., Phadnis N.V., Suryanarayanan R., *Powder Diffraction*, 2001, 16(1), 20.

MS20 P11

Structural investigation of new insulin derivative at room temperature Biserka Prugovečki^a, Stjepan Prugovečki^b, Detlef Beckers^b, Dubravka Matković-Čalogović^a ^aDepartment of Chemistry, University of Zagreb, Croatia. ^bPanalytical B.V., Almelo, The Netherlands. E-mail: biserka@chem.pmf.hr

Keywords: insulin, protein crystallography, powder diffraction

Insulin is a hormone protein that regulates carbohydrate metabolism and it also takes part in the metabolism of fat and proteins. It is used medically in patients with Type 1 diabetes mellitus. Occasionally some patients with Type 2 diabetes mellitus also require insulin.

Owing to its crucial metabolic role and its pharmaceutical importance many structural studies on chemically and genetically modified insulins have been done.

We will present the results of our investigation on human bromo-derivative of insulin. Both single crystal and powder diffraction data were collected on laboratory instruments at room temperature. The investigated insulin derivative belongs to the T₃R₃^f rhombohedral form [1] with cell parameters $a = 80.96 \text{ \AA}$ and $c = 37.30 \text{ \AA}$. The unit cell parameter c at room temperature is two times smaller in comparison to the one at 100 K [2]. Coordination of zinc ions and conformation of insulin molecule will be discussed.

[1] Kaarsholm, N. C., Ko H. C., Dunn M. F. *Biochemistry*, 1989, 28, 4427.

[2] I. Đilović I. et al., unpublished results.

MS20 P12

Structure solution and metastable zone width experiments of a tri-substituted aromatic compound. Andrew O'Neill, Chick C. Wilson, WestCHEM, Department of Chemistry, University of Glasgow, Glasgow G12 8QQ, Alastair J. Florence, Department of Pharmaceutical Sciences, University of Strathclyde, 27 Taylor Street, Glasgow G4 0NR.

Keywords: polymorphism, crystal nucleation, metastable zone

Polymorphism in molecular crystals is the ability of a substance to exist in different crystal packing arrangements [1]. Understanding the phenomenon of polymorphism has become an increasingly important challenge, particularly in the pharmaceutical industry, where there would be considerable advantages were it possible to be able to identify which compounds will be likely to display different polymorphic forms from a knowledge of molecular structure alone. We are part of the UK Research Councils' Basic Technology CPOSS project (Control and Prediction of the Organic Solid State), which has been set up to tackle this problem.

Whilst accurate thermodynamic models are available to predict polymorphism, they do not currently accommodate kinetic factors such as nucleation and crystal growth resulting in an inclination to overestimate the tendency to polymorphism [2]. Nucleation is the initial process leading to the growth of crystals. Due to rapid onset, nucleation studies have proven to be a highly challenging area to study experimentally. However, an understanding of the process is essential as increasing numbers of important materials are found to exhibit polymorphism, even when grown under seemingly identical conditions.

This poster will describe the experimental approach to study model systems, the identification and structure solution of suitable systems and establishment of conditions that are suitable for subsequent examination by scattering techniques.

Methyl 2,5-dibromobenzoate (C₈H₆O₂ Br₂), a tri-substituted aromatic compound, with a previously unsolved crystal structure, has been selected for the investigations. The choice of this material is governed by